


Instrumentation available in the Microscopy Core Facility

Observation of samples in Petri dishes or multi-well plates



Axiovert S100TV Inverted Microscope

The Axiovert S100TV inverted microscope from ZEISS is equipped with a motorized MS-2000 controller board from ASI Imaging (Applied Scientific Instrumentation). Different racks are available to deposit kneaders of various sizes or multi-well plates. The Retiga EXi camera is used to acquire monochromatic images of samples. Several manual filter cubes are available.



DMIRE Inverted Microscope

The LEICA DMIRE microscope allows the observation of different samples in transmitted light or by fluorescence. Several manual filter cubes are available. If required, images can be acquired by moving a Retiga EXi camera from a nearby system.